



AN 373

Process Control Using SIMS for Compound Semiconductors

May 7, 2007 (Version 3.0)

Discussion

The measurement of dopant and impurity concentrations as a function of depth in compound semiconductors is often accomplished using Secondary Ion Mass Spectrometry (SIMS), because of its very low detection limits and excellent depth resolution. SIMS is a key analytical technique for research and development of compound semiconductor epitaxy and devices, but does it have the capability for process control?

Process control, or monitoring, typically requires a tight precision over a long period of time. The required measurement precision is dictated by the required process tolerance, i.e., the measurement precision must obviously be tighter (or smaller) than the process tolerance. The ratio of the process tolerance (PT), in terms of a standard deviation, to the measurement precision (MP) is sometimes called the PT/MP ratio, and depending upon one's use of statistical control in manufacturing the desired PT/MP ratio varies from 3x to 10x. (For example, the International Technology Roadmap for Semiconductors, ITRS, targets metrology for a ratio of 3x minimum to a desired 10x if possible. However, there are some metrologies in use in the silicon semiconductor industry that do not meet even the 3x requirement.)

An indication of what is possible with SIMS for compound semiconductors is revealed in the following example. Figure 1 shows a SIMS depth profile of Si, C, O, H, In, Ga, and Al in a GaAs HBT wafer. This is a planar structure, i.e., not a device, so many "samples" can be taken from this wafer over time for SIMS measurements. This is important since SIMS is a destructive technique. The assumption is that the HBT structure is uniform across the wafer.

The important point is that this sample, a CONTROL SAMPLE, is measured after a reference material is analyzed to determine relative sensitivity factors. The CONTROL SAMPLE has been measured many times over a period of years, a longer period of time typically than the life of any one compound semiconductor process. The precision determined from these measurements include variations due to the measurement process plus any non-uniformity across the wafer.

Figures 2, 3 and 4 are long term charts of different parameters taken from the SIMS measurement of the CONTROL SAMPLE. For example, Figure 2 shows the Carbon Base Doping, the Silicon Top Contact Doping, and the Silicon Sub-Collector Doping versus time. The relative standard deviation (RSD), as defined by the ratio of the one standard deviation to the average value, is shown in the chart for these parameters. The RSD varies from 4.7% to 6.6%. Also shown in Figure 2 is the precision for measurement of the Al atom fraction. Even for a key matrix component, the RSD is 2.9%. Figure 3 shows the data for the Silicon Collector




Figure 1. Control Sample for Long-Term Precision Check




Figure 2. Long-Term Precision Statistics

Doping, and this has an RSD of 8.2%. Note that the higher RSD corresponds to a much lower doping level (4×10^{16} atoms/cm³) than the parameters in Figure 2 ($>10^{18}$ atoms/cm³). This illustrates a key point that the SIMS precision varies with concentration level.

Because of the excellent depth resolution of SIMS profiles it is also possible to accurately monitor layer thickness or junction depth. This is shown in Figure 4 for the epi thickness and Base/Collector Junction Depth. The RSDs of these parameters are 1.2% and 1.5% respectively.

In summary, SIMS can be used for long-term process control or monitoring of compound semiconductors. Whether SIMS is the measurement of choice for a particular application will depend upon a determination of the SIMS precision and upon the process tolerance of the manufacturing process.




Figure 3. Long-Term Precision Statistics




Figure 4. Long-Term Precision Statistics

United States Locations

Tempe, Arizona
+1 480 239 0602 info.az@eaglabs.com
+1 602 470 2655 fax

Sunnyvale, California
810 Kifer Road
+1 408 530 3500 info.ca@eaglabs.com
+1 408 530 3501 fax

1135 E Arques Avenue
+1 408 738 3033
+1 408 738 3035 fax

785 Lucerne Drive
+1 408 737 3892
+1 408 737 3916 fax

Peabody, Massachusetts
+1 978 278 9500 info.ma@eaglabs.com
+1 978 278 9501 fax

Chanhassen, Minnesota
+1 952 828 6411 info.mn@eaglabs.com
+1 952 828 6449 fax

East Windsor, New Jersey
+1 609 371 4800 info.nj@eaglabs.com
+1 609 371 5666 fax

Syracuse, New York
+1 315 431 9900 info.ny@eaglabs.com
+1 315 431 9800 fax

Raleigh, North Carolina
+1 919 829 7041 info.nc@eaglabs.com
+1 919 829 5518 fax

Round Rock, Texas
+1 512 671 9500 info.tx@eaglabs.com
+1 512 671 9501 fax

International Locations

Shanghai, China
+ 86 21 6879 6088 info.cn@eaglabs.com
+ 86 21 6879 9086 fax

Tournefeuille, France
+ 33 5 61 73 15 29 info.fr@eaglabs.com
+ 33 5 61 73 15 67 fax

Frankfurt, Germany
+ 49 (0) 693053213 info.de@eaglabs.com
+ 49 (0) 69307941 fax

Tokyo, Japan
+ 81 3 5396 0531 info.jp@eaglabs.com
+ 81 3 5396 1930 fax

HsinChu, Taiwan
+ 886 3 5632303 info.tw@eaglabs.com
+ 886 3 5632306 fax

Uxbridge, United Kingdom
+ 44 (0) 1895 811194 info.uk@eaglabs.com
+ 44 (0) 1895 810350 fax