

# Analytical Solutions

for BioTechnology

May 9, 2007 (Version 2.0)

## BN 1504

### Measuring the Morphological Structure of Porous Materials

#### Introduction

Porous materials are important in many industries including pharmaceuticals, tissue engineering, and biomedical technology. Applications for these materials include protein separation, drug delivery, surface modification, scaffolding for tissue growth, and support material in implant technology. The ability to determine parameters such as pore size, pore shape, and interconnectivity are critical in research and development in the areas of filtration, permeability, tissue engineering, cell sustainability, and timed release of active pharmaceutical materials.

#### Discussion

Current technology produces porous structures with pore sizes and fiber diameters on the order of nanometers. The low-voltage field-emission scanning electron microscope (FE SEM) is ideally suited for characterizing the morphology of porous materials. FE SEM provides the combination of high spatial resolution and depth of field which are necessary for imaging these materials. In addition the FE-SEM can be operated at low accelerating voltages without a substantial loss of resolution. Operating at 1keV minimizes charging, beam damage, and the necessity for additional sample preparation. Imaging at lower voltages also provides more sensitivity to the delicate surface features common in these systems. FE-SEM images of the top and bottom exterior surfaces and the cross-sectioned interior pore structures of a filter membrane are shown in Figure 1.





Figure 1. FE-SEM images of (A) top exterior surface, (B) interior cross-section, and (C) bottom exterior surface of a filter membrane.



---

#### United States Locations

Tempe, Arizona  
+1 480 239 0602 [info.az@eaglabs.com](mailto:info.az@eaglabs.com)  
+1 602 470 2655 fax

Sunnyvale, California  
810 Kifer Road  
+1 408 530 3500 [info.ca@eaglabs.com](mailto:info.ca@eaglabs.com)  
+1 408 530 3501 fax

1135 E Arques Avenue  
+1 408 738 3033  
+1 408 738 3035 fax

785 Lucerne Drive  
+1 408 737 3892  
+1 408 737 3916 fax

Peabody, Massachusetts  
+1 978 278 9500 [info.ma@eaglabs.com](mailto:info.ma@eaglabs.com)  
+1 978 278 9501 fax

Chanhassen, Minnesota  
+1 952 828 6411 [info.mn@eaglabs.com](mailto:info.mn@eaglabs.com)  
+1 952 828 6449 fax

East Windsor, New Jersey  
+1 609 371 4800 [info.nj@eaglabs.com](mailto:info.nj@eaglabs.com)  
+1 609 371 5666 fax

Syracuse, New York  
+1 315 431 9900 [info.ny@eaglabs.com](mailto:info.ny@eaglabs.com)  
+1 315 431 9800 fax

Raleigh, North Carolina  
+1 919 829 7041 [info.nc@eaglabs.com](mailto:info.nc@eaglabs.com)  
+1 919 829 5518 fax

Round Rock, Texas  
+1 512 671 9500 [info.tx@eaglabs.com](mailto:info.tx@eaglabs.com)  
+1 512 671 9501 fax

#### International Locations

Shanghai, China  
+ 86 21 6879 6088 [info.cn@eaglabs.com](mailto:info.cn@eaglabs.com)  
+ 86 21 6879 9086 fax

Tournefeuille, France  
+ 33 5 61 73 15 29 [info.fr@eaglabs.com](mailto:info.fr@eaglabs.com)  
+ 33 5 61 73 15 67 fax

Frankfurt, Germany  
+ 49 (0) 693053213 [info.de@eaglabs.com](mailto:info.de@eaglabs.com)  
+ 49 (0) 69307941 fax

Tokyo, Japan  
+ 81 3 5396 0531 [info.jp@eaglabs.com](mailto:info.jp@eaglabs.com)  
+ 81 3 5396 1930 fax

HsinChu, Taiwan  
+ 886 3 5632303 [info.tw@eaglabs.com](mailto:info.tw@eaglabs.com)  
+ 886 3 5632306 fax

Uxbridge, United Kingdom  
+ 44 (0) 1895 811194 [info.uk@eaglabs.com](mailto:info.uk@eaglabs.com)  
+ 44 (0) 1895 810350 fax